

Anti-TBXA2R Antibody (5B707)

Product Details

Ig Type:	IgG
Reactivity:	Human,Mouse,Rat
Conjugation:	Unconjugated
Clone:	5B707
Purification:	Affinity-chromatography

Applications

Verified Activity:	Western blot analysis of Thromboxane A2 receptor expression in human platelet cell lysate.
Application:	WB
Recommended	WB: 1:1000-1:2000

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	A synthesized peptide: human TBXA2R
Synonyms:	TXA2R;BDPLT13;TP receptor;Tbxa2r

Research Background

Receptor for thromboxane A2 (TXA2), a potent stimulator of platelet aggregation.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481