

Anti-Thyroglobulin Antibody (6Q665)

Product Details

Ig Type:	IgG
Reactivity:	Human
Conjugation:	Unconjugated
Clone:	6Q665
Purification:	Affinity-chromatography

Applications

Verified Activity:	Western blot analysis of Thyroglobulin expression in Human thyroid lysate.
Application:	IHC,IP,WB
Recommended	IHC: 1:100-1:200; IP: 1:20-1:50; WB: 1:1000-1:2000

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	A synthesized peptide: human Thyroglobulin
Synonyms:	TGN;AITD3;Tg;hTG;TDH3

Research Background

Acts as a substrate for the production of iodinated thyroid hormones thyroxine (T4) and triiodothyronine (T3)

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481