

Anti-Trypsin 2 Antibody (5C6)

Product Details

Ig Type:	IgG
Reactivity:	Human,Mouse,Rat
Conjugation:	Unconjugated
Clone:	5C6
Purification:	Affinity-chromatography

Applications

Verified Activity:	Western blot analysis of PRSS2 expression in Human pancreas lysate.
Application:	WB
Recommended	WB: 1:1000-1:2000

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	A synthesized peptide: human Trypsin 2
Synonyms:	Protease;Pancreatic anionic trypsinogen;serine 2;PRSS2;TRYP8;Trypsin2

Research Background

Part of the insoluble cornified cell envelope (CE) of stratified squamous epithelia.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481