

Anti-TSH beta Antibody (4017)

Product Details

Ig Type:	IgG
Reactivity:	Human
Conjugation:	Unconjugated
Clone:	4017
Purification:	Affinity-chromatography

Applications

Verified Activity:	Western blot analysis of TSH beta expression in Human pituitary lysate.
Application:	IHC,WB
Recommended	IHC: 1:100-1:200; WB: 1:1000-1:2000

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	A synthesized peptide: human TSH beta
Synonyms:	CHNG4;TSH β ;TSH b;TSHB;TSH β

Research Background

Indispensable for the control of thyroid structure and metabolism.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481