

Anti-TXNDC12 Antibody (5W255)

Product Details

Ig Type:	IgG
Reactivity:	Human,Mouse,Rat
Conjugation:	Unconjugated
Clone:	5W255
Purification:	Affinity-chromatography

Applications

Verified Activity:	Western blot analysis of ERp19 expression in 293T cell lysate.
Application:	ICC/IF,IHC,IHC-Fr,WB
Recommended	ICC/IF: 1:50-1:200; IHC: 1:100-1:200; IHC-Fr: 1:50-1:200; WB: 1:1000-1:2000

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	A synthesized peptide: human TXNDC12
Synonyms:	ERP19;hAG1;ER protein 19;ERP16;hTLP19;PDIA16;AG1;AGR1;TLP19;TXNDC 12;ER protein 18;ERP 18

Research Background

Protein-disulfide reductase of the endoplasmic reticulum that promotes disulfide bond formation in client proteins through its thiol-disulfide oxidase activity.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481