

Anti-UBA7 Antibody (3C20)

Product Details

Ig Type:	IgG
Reactivity:	Human,Mouse,Rat
Conjugation:	Unconjugated
Clone:	3C20
Purification:	Affinity-chromatography

Applications

Verified Activity:	Western blot analysis of Ube1L / UBA7 in HepG2 cell lysate.
Application:	FCM,IHC,WB
Recommended	FCM: 1:20-1:100; IHC: 1:100-1:200; WB: 1:1000-1:2000

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	A synthesized peptide: human UBA7
Synonyms:	UBA 7;Ubiquitin activating enzyme 2;UBE2;UBE1L;Ubiquitin-activating enzyme 7;UBA1B; Ubiquitin-activating enzyme E1 homolog

Research Background

E1-activating enzyme that catalyzes the covalent conjugation of the ubiquitin-like protein product of ISG15 to additional interferon stimulated proteins (ISGs) as well as other cellular proteins such as P53 in a process termed protein ISGylation.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481