

Anti-UBE4A Antibody (6C483)

Product Details

| | |
|---------------|-------------------------|
| Ig Type: | IgG |
| Reactivity: | Human |
| Conjugation: | Unconjugated |
| Clone: | 6C483 |
| Purification: | Affinity-chromatography |

Applications

| | |
|--------------------|---|
| Verified Activity: | Western blot analysis of UBE4A expression in 293 cell lysate. |
| Application: | FCM,ICC/IF,WB |
| Recommended | FCM: 1:20-1:100; ICC/IF: 1:50-1:200; WB: 1:1000-1:2000 |

Properties

| | |
|----------------------|---|
| Stability & Storage: | Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles. |
| Shipping: | Shipping with blue ice. |

Antigen Details

| | |
|------------|------------------------------------|
| Immunogen: | A synthesized peptide: human UBE4A |
| Synonyms: | UFD2;UBE 4A;UBOX2;UFD2b |

Research Background

Ubiquitin-protein ligase that probably functions as an E3 ligase in conjunction with specific E1 and E2 ligases. May also function as an E4 ligase mediating the assembly of polyubiquitin chains on substrates ubiquitinated by another E3 ubiquitin ligase. Mediates 'Lys-48'-linked polyubiquitination of substrates.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481