

## Anti-VAV1/2/3 Antibody (4M772)

### Product Details

Ig Type:	IgG
Reactivity:	Human
Conjugation:	Unconjugated
Clone:	4M772
Purification:	Affinity-chromatography

### Applications

Verified Activity:	Western blot analysis of Vav1/2/3 expression in Jurkat cell lysate.
Application:	FCM,ICC/IF,IP,WB
Recommended	FCM: 1:20-1:100; ICC/IF: 1:50-1:200; IP: 1:20-1:50; WB: 1:1000-1:2000

### Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

### Antigen Details

Immunogen:	A synthesized peptide: human VAV1
Synonyms:	Protein vav 2;p95;vav-T;VAV;Protein vav 3;Protein vav 1;VAV 1/2/3;p95vav

### Research Background

Couples tyrosine kinase signals with the activation of the Rho/Rac GTPases, thus leading to cell differentiation and/or proliferation.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481