

## Anti-VSX2 Antibody (1Q587)

### Product Details

Ig Type:	IgG
Reactivity:	Human,Mouse,Rat
Conjugation:	Unconjugated
Clone:	1Q587
Purification:	Affinity-chromatography

### Applications

Verified Activity:	Western blot analysis of CHX10 expression in Y79 cell lysate.
Application:	FCM,WB
Recommended	FCM: 1:20-1:100; WB: 1:1000-1:2000

### Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

### Antigen Details

Immunogen:	A synthesized peptide: human VSX2
Synonyms:	CHX 10;Vsx2;VSX 2;HOX10;RET1;MCOPCB3;MCOP2

### Research Background

Plays a significant role in the specification and morphogenesis of the sensory retina. May also participate in the development of the cells of the inner nuclear layer, particularly bipolar cells.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481