

Anti-WAS Antibody (6B331)

Product Details

Ig Type:	IgG
Reactivity:	Human,Mouse,Rat
Conjugation:	Unconjugated
Clone:	6B331
Purification:	Affinity-chromatography

Applications

Verified Activity:	Western blot analysis of WASP expression in Jurkat cell lysate.
Application:	FCM,ICC/IF,IP,WB
Recommended	FCM: 1:20-1:100; ICC/IF: 1:50-1:200; IP: 1:20-1:50; WB: 1:1000-1:2000

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	A synthesized peptide: human WAS
Synonyms:	Was;THC1;THC;SCNX;WASp;U42471;IMD2

Research Background

Effector protein for Rho-type GTPases, providing a link with the Arp2/3 complex that regulates the structure and dynamics of the actin cytoskeleton. Important for efficient actin polymerization. Possible regulator of lymphocyte and platelet function.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481