

## Anti-Human/Mouse CD95 Antibody (HFE7A)

## Chemical Properties

CAS No. :

Formula:

Molecular Weight: 150 kDa

Store at low temperature

Storage: Store at -20°C

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	Anti-Human/Mouse CD95 Antibody (HFE7A) is a monoclonal inhibitory antibody that recognizes both human and mouse CD95 (Fas) and induces apoptosis in thymocytes.
Targets(IC50)	Immunology/Inflammation related
In vivo	Anti-Human/Mouse CD95 Antibody (HFE7A) can induce apoptosis in mouse thymocytes.

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.0067 mL	0.0333 mL	0.0667 mL
5 mM	0.0013 mL	0.0067 mL	0.0133 mL
10 mM	0.0007 mL	0.0033 mL	0.0067 mL
50 mM	0.0001 mL	0.0007 mL	0.0013 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481