

4",5"-dehydroisopsoralidin

Chemical Properties


CAS No. : 65639-51-2

Formula: C₂₀H₁₄O₅

Molecular Weight: 334.32

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	4",5"-dehydroisopsoralidin is a natural product.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.9911 mL	14.9557 mL	29.9115 mL
5 mM	0.5982 mL	2.9911 mL	5.9823 mL
10 mM	0.2991 mL	1.4956 mL	2.9911 mL
50 mM	0.0598 mL	0.2991 mL	0.5982 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481