


## 2,4-Dibromoaniline

## Chemical Properties

CAS No. :	615-57-6
Formula:	C <sub>6</sub> H <sub>5</sub> Br <sub>2</sub> N
Molecular Weight:	250.919
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	2, 4-dibromoaniline belongs to the 2-bromoaniline class of organic compounds. 2, 4-dibromoaniline is often used as a model substance for studying the chemical property changes of 2-bromoaniline compounds after the introduction of polyhalogen atoms.
-------------	--

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.9853 mL	19.9267 mL	39.8533 mL
5 mM	0.7971 mL	3.9853 mL	7.9707 mL
10 mM	0.3985 mL	1.9927 mL	3.9853 mL
50 mM	0.0797 mL	0.3985 mL	0.7971 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481