


CFMTI

Chemical Properties

CAS No. :	864864-17-5
Formula:	C ₁₉ H ₁₆ FN ₅ O
Molecular Weight:	349.36
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	CFMTI is a potent and selective metabotropic glutamate receptor (mGluR) 1 allosteric antagonist (IC ₅₀ : 2.6 nM). The selectivity of CFMTI to mGluR1 is higher 2000-fold than mGluR5.
Targets(IC ₅₀)	GluR

Solubility Information

Solubility	DMSO: 4.62 mg/mL (13.22 mM), Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.8624 mL	14.3119 mL	28.6238 mL
5 mM	0.5725 mL	2.8624 mL	5.7248 mL
10 mM	0.2862 mL	1.4312 mL	2.8624 mL
50 mM	0.0572 mL	0.2862 mL	0.5725 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Satow A et al. Unique antipsychotic activities of the selective metabotropic glutamate receptor 1 allosteric antagonist 2-cyclopropyl-5-[1-(2-fluoro-3-pyridinyl)-5-methyl-1H-1,2,3-triazol-4-yl]-2,3-dihydro-1H-isoindol-1-one. J Pharmacol Exp Ther. 2009 Jul;330(1):179-90.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481