

## (2R,3S)-Pteroside C

## Chemical Properties

CAS No. :	68399-16-6
Formula:	C <sub>20</sub> H <sub>28</sub> O <sub>8</sub>
Molecular Weight:	396.43
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

## Biological Description

Description	(2R,3S)-Pteroside C is a useful organic compound for research related to life sciences and the catalog number is T124409.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.5225 mL	12.6126 mL	25.2251 mL
5 mM	0.5045 mL	2.5225 mL	5.045 mL
10 mM	0.2523 mL	1.2613 mL	2.5225 mL
50 mM	0.0505 mL	0.2523 mL	0.5045 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481