


## 4-(Phenyldiazenyl)benzoic acid

## Chemical Properties

CAS No. :	1562-93-2
Formula:	C <sub>13</sub> H <sub>10</sub> N <sub>2</sub> O <sub>2</sub>
Molecular Weight:	226.23
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	4-(Phenyldiazenyl)benzoic acid is a photosensitive and photoswitchable TRPA1 agonist used as a pharmacological tool to study pain signaling.
Targets(IC50)	Others,TRP/TRPV Channel

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.4203 mL	22.1014 mL	44.2028 mL
5 mM	0.8841 mL	4.4203 mL	8.8406 mL
10 mM	0.442 mL	2.2101 mL	4.4203 mL
50 mM	0.0884 mL	0.442 mL	0.8841 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481