

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.2751 mL	11.3755 mL	22.7511 mL
5 mM	0.455 mL	2.2751 mL	4.5502 mL
10 mM	0.2275 mL	1.1376 mL	2.2751 mL
50 mM	0.0455 mL	0.2275 mL	0.455 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Christopher D Kane, et al. Molecular characterization of novel and selective peroxisome proliferator-activated receptor alpha agonists with robust hypolipidemic activity in vivo. Mol Pharmacol. 2009 Feb;75(2):296-306.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481