


## Daturaturin A

## Chemical Properties

CAS No. :	133360-51-7
Formula:	C34H48O10
Molecular Weight:	616.748
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Daturaturin A is a natural product of Datura, Solanaceae. The catalog number is TN3776 and the CAS number is 133360-51-7. Daturaturin A can be used as a reference standard.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.6214 mL	8.107 mL	16.214 mL
5 mM	0.3243 mL	1.6214 mL	3.2428 mL
10 mM	0.1621 mL	0.8107 mL	1.6214 mL
50 mM	0.0324 mL	0.1621 mL	0.3243 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481