

SJB2-043

## Chemical Properties


CAS No. : 63388-44-3

Formula: C<sub>17</sub>H<sub>9</sub>NO<sub>3</sub>

Molecular Weight: 275.26

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	SJB2-043, an inhibitor of USP1-UAF1, has an IC <sub>50</sub> of 544 nM.
Targets(IC <sub>50</sub> )	DUB

## Solubility Information

Solubility	DMSO: 3.33 mg/mL (12.1 mM), Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	---

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.6329 mL	18.1646 mL	36.3293 mL
5 mM	0.7266 mL	3.6329 mL	7.2659 mL
10 mM	0.3633 mL	1.8165 mL	3.6329 mL
50 mM	0.0727 mL	0.3633 mL	0.7266 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

## Reference

Jussara, Maria, Goncalves, et al. The Role of the Complex USP1/WDR48 in Differentiation and Proliferation Processes in Cancer Stem Cells.[J]. Current stem cell research & therapy, 2017.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481