

Gentianose

Chemical Properties

CAS No. : 25954-44-3


Formula: C₁₈H₃₂O₁₆

Molecular Weight: 504.44

Keep away from moisture, Keep away from direct sunlight

Storage: Pure form: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	Gentianose is a carbohydrate reserve derived from perennial yellow gentian storage roots.
Targets(IC50)	Others

Solubility Information

Solubility	DMSO: 30 mg/mL (59.47 mM), Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.9824 mL	9.912 mL	19.824 mL
5 mM	0.3965 mL	1.9824 mL	3.9648 mL
10 mM	0.1982 mL	0.9912 mL	1.9824 mL
50 mM	0.0396 mL	0.1982 mL	0.3965 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286

E_mail: info@targetmol.com

Address: 34 Washington Street, Wellesley Hills, MA 02481