

Azido-PEG7-alcohol

Chemical Properties

CAS No. : 1274892-60-2


Formula: C14H29N3O7

Molecular Weight: 351.4

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	Azido-PEG7-alcohol is a PEG-based linker for PROTACs that connects two essential ligands, facilitating the formation of PROTAC molecules and enabling selective protein degradation via the ubiquitin-proteasome system within cells (brackets).
Targets(IC50)	Others,PROTAC Linker
In vitro	PROTACs comprise two distinct ligands joined by a linker: one ligand binds to an E3 ubiquitin ligase, and the other to the target protein. They leverage the intracellular ubiquitin-proteasome system to selectively degrade target proteins [1].

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.8458 mL	14.2288 mL	28.4576 mL
5 mM	0.5692 mL	2.8458 mL	5.6915 mL
10 mM	0.2846 mL	1.4229 mL	2.8458 mL
50 mM	0.0569 mL	0.2846 mL	0.5692 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

An S, et al. Small-molecule PROTACs: An emerging and promising approach for the development of targeted therapy drugs. EBioMedicine. 2018 Oct;36:553-562

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481