

$\alpha$ -CGRP (rat) acetate

## Chemical Properties

CAS No. :

Formula:

Molecular Weight: 3806.30

Store at low temperature,Keep away from moisture

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	$\alpha$ -CGRP (rat) acetate is a calcitonin gene-related peptide (CGRP). $\alpha$ -CGRP (rat) acetate is a potent vasodilator.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.2627 mL	1.3136 mL	2.6272 mL
5 mM	0.0525 mL	0.2627 mL	0.5254 mL
10 mM	0.0263 mL	0.1314 mL	0.2627 mL
50 mM	0.0053 mL	0.0263 mL	0.0525 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

## Reference

S D Brain, et al. Calcitonin gene-related peptide is a potent vasodilator. Nature. 1985 Jan 3-9;313(5997):54-6.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481