

Zeylasteral

Chemical Properties


CAS No. :

Formula: C₃₀H₃₈O₆

Molecular Weight: 494.628

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	Zeylasteral is a natural product that can be used as a reference standard.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.0217 mL	10.1086 mL	20.2171 mL
5 mM	0.4043 mL	2.0217 mL	4.0434 mL
10 mM	0.2022 mL	1.0109 mL	2.0217 mL
50 mM	0.0404 mL	0.2022 mL	0.4043 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481