


S-2-(4-Aminobenzyl)-diethylenetriamine pentaacetic acid TFA salt, p-NH₂-Bn-DTPA-TFA

Chemical Properties

CAS No. : 102650-29-3
 Formula: C₂₁H₃₀N₄O₁₀
 Molecular Weight: 498.48
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	Super-chelating agent & MRI contrast agent
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.0061 mL	10.0305 mL	20.061 mL
5 mM	0.4012 mL	2.0061 mL	4.0122 mL
10 mM	0.2006 mL	1.003 mL	2.0061 mL
50 mM	0.0401 mL	0.2006 mL	0.4012 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481