


Ganoderic acid SZ

Chemical Properties

CAS No. :	865543-37-9
Formula:	C ₃₀ H ₄₄ O ₃
Molecular Weight:	452.679
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Ganoderic acid SZ is a natural product. The catalog number is TN4104 and the CAS number is 865543-37-9. Ganoderic acid SZ can be used as a reference standard.
Targets(IC50)	Others,HMG-CoA Reductase

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.2091 mL	11.0453 mL	22.0907 mL
5 mM	0.4418 mL	2.2091 mL	4.4181 mL
10 mM	0.2209 mL	1.1045 mL	2.2091 mL
50 mM	0.0442 mL	0.2209 mL	0.4418 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481