

SR7826

Chemical Properties


CAS No. : 1219728-20-7

Formula: C₂₂H₂₁N₅O₂

Molecular Weight: 387.43

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	SR7826 is a selective LIMK inhibitor.
Targets(IC50)	Others,LIM Kinase

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.5811 mL	12.9056 mL	25.8111 mL
5 mM	0.5162 mL	2.5811 mL	5.1622 mL
10 mM	0.2581 mL	1.2906 mL	2.5811 mL
50 mM	0.0516 mL	0.2581 mL	0.5162 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Yu Q, Gratzke C, Wang Y, Herlemann A, Sterr CM, Rutz B, Ciotkowska A, Wang X, Strittmatter F, Stief CG, Hennenberg M. Inhibition of human prostate smooth muscle contraction by the LIM kinase inhibitors, SR7826 and LIMKi3. Br J Pharmacol. 2018 Jun;175(11):2077-2096. doi: 10.1111/bph.14201. Epub 2018 Apr 29. PubMed PMID: 29574791; PubMed Central PMCID: PMC5978953.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481