

Cepharadione B

Chemical Properties

CAS No. :	55610-02-1
Formula:	C ₁₉ H ₁₅ NO ₄
Molecular Weight:	321.33
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

Biological Description

Description	Cepharadione B is a useful organic compound for research related to life sciences and the catalog number is T124661.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.1121 mL	15.5603 mL	31.1207 mL
5 mM	0.6224 mL	3.1121 mL	6.2241 mL
10 mM	0.3112 mL	1.556 mL	3.1121 mL
50 mM	0.0622 mL	0.3112 mL	0.6224 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481