

Compound N038-0206

Chemical Properties


CAS No. :

Formula: C₁₈H₂₃NO₄

Molecular Weight: 317.39

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	Compound N038-0206 is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.1507 mL	15.7535 mL	31.507 mL
5 mM	0.6301 mL	3.1507 mL	6.3014 mL
10 mM	0.3151 mL	1.5753 mL	3.1507 mL
50 mM	0.063 mL	0.3151 mL	0.6301 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481