

(4-Oxo-4H-quinazolin-3-yl)-acetic acid

Chemical Properties

CAS No. : 14663-53-7


Formula: C₁₀H₈N₂O₃

Molecular Weight: 204.1821

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	(4-Oxo-4H-quinazolin-3-yl)-acetic acid is an alkyl chain-derived PROTAC linker applicable for the synthesis of PROTACs [1].
Targets(IC50)	Others,PROTAC Linker
In vitro	PROTACs, consisting of two distinct ligands joined by a linker—one binding to an E3 ubiquitin ligase and the other to the target protein—utilize the intracellular ubiquitin-proteasome system to selectively degrade target proteins[1].

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.8976 mL	24.4882 mL	48.9764 mL
5 mM	0.9795 mL	4.8976 mL	9.7953 mL
10 mM	0.4898 mL	2.4488 mL	4.8976 mL
50 mM	0.098 mL	0.4898 mL	0.9795 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

An S, et al. Small-molecule PROTACs: An emerging and promising approach for the development of targeted therapy drugs. EBioMedicine. 2018 Oct;36:553-562.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481