

## Stachysterone D

## Chemical Properties


CAS No. :

Formula: C<sub>27</sub>H<sub>42</sub>O<sub>6</sub>

Molecular Weight: 462.627

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Stachysterone D is a useful organic compound for research related to life sciences and the catalog number is T125165.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.1616 mL	10.8078 mL	21.6155 mL
5 mM	0.4323 mL	2.1616 mL	4.3231 mL
10 mM	0.2162 mL	1.0808 mL	2.1616 mL
50 mM	0.0432 mL	0.2162 mL	0.4323 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481