

HbF inducer-1

Chemical Properties

CAS No. :

Formula: C₁₈H₁₉N₃O₃

Molecular Weight: 325.36

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

| | |
|---------------|---|
| Description | HbF inducer-1 is a fetal hemoglobin inducer which is orally bioavailable. |
| Targets(IC50) | Others,Histone Methyltransferase |

Preparing Stock Solutions

| | 1mg | 5mg | 10mg |
|-------|-----------|------------|------------|
| 1 mM | 3.0735 mL | 15.3676 mL | 30.7352 mL |
| 5 mM | 0.6147 mL | 3.0735 mL | 6.147 mL |
| 10 mM | 0.3074 mL | 1.5368 mL | 3.0735 mL |
| 50 mM | 0.0615 mL | 0.3074 mL | 0.6147 mL |

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481