

## Compound T130720

## Chemical Properties


CAS No. :

Formula: C<sub>29</sub>H<sub>39</sub>NO<sub>11</sub>S

Molecular Weight: 609.69

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	2-acetamido-3-[(8,9,14,20-tetrahydroxy-3,7,13,16-tetramethyl-2,19-dioxo-18,21-dioxaheptacyclo[11.6.2.0 <sup>1,14</sup> .0 <sup>3,12</sup> .0 <sup>4,9</sup> .0 <sup>17,20</sup> ]]henicos-10-en-17-yl)methyl)sulfanyl]propanoic acid is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	---

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.6402 mL	8.2009 mL	16.4018 mL
5 mM	0.328 mL	1.6402 mL	3.2804 mL
10 mM	0.164 mL	0.8201 mL	1.6402 mL
50 mM	0.0328 mL	0.164 mL	0.328 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481