

Sanggenol P

Chemical Properties


CAS No. : 1351931-30-0

Formula: C₃₀H₃₆O₆

Molecular Weight: 492.6

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	Sanggenol P is a natural product of Morus, Moraceae. The catalog number is TN4940 and the CAS number is 1351931-30-0. Sanggenol P can be used as a reference
Targets(IC50)	Others,HBV

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.030 mL	10.1502 mL	20.3004 mL
5 mM	0.406 mL	2.030 mL	4.0601 mL
10 mM	0.203 mL	1.015 mL	2.030 mL
50 mM	0.0406 mL	0.203 mL	0.406 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481