


C 175

## Chemical Properties

CAS No. : 6237-31-6  
 Formula: C<sub>19</sub>H<sub>24</sub>Cl<sub>2</sub>N<sub>4</sub>O<sub>2</sub>  
 Molecular Weight: 411.33  
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
 Actual storage temperature shall be subject to the COA.



## Biological Description

Description	C 175 is a bioactive agent.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.4311 mL	12.1557 mL	24.3114 mL
5 mM	0.4862 mL	2.4311 mL	4.8623 mL
10 mM	0.2431 mL	1.2156 mL	2.4311 mL
50 mM	0.0486 mL	0.2431 mL	0.4862 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481