

Compound N005-0008

## Chemical Properties


CAS No. :

Formula: C<sub>15</sub>H<sub>15</sub>N<sub>3</sub>O<sub>3</sub>S

Molecular Weight: 317.36

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

|               |  |
|---------------|--|
| Description   | Compound N005-0008 is a useful organic compound for research related to life sciences and the catalog number is T131377. |
| Targets(IC50) | Others   |

## Preparing Stock Solutions

|       | 1mg       | 5mg       | 10mg      |
|-------|-----------|-----------|-----------|
| 1 mM  | 3.151 mL  | 15.755 mL | 31.510 mL |
| 5 mM  | 0.6302 mL | 3.151 mL  | 6.302 mL  |
| 10 mM | 0.3151 mL | 1.5755 mL | 3.151 mL  |
| 50 mM | 0.063 mL  | 0.3151 mL | 0.6302 mL |

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481