


2'-O-Methyl-2,5'-anhydrouridine

Chemical Properties

CAS No. : 175471-64-4
 Formula: C₁₀H₁₂N₂O₅
 Molecular Weight: 240.21
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	2'-O-Methyl-2,5'-anhydrouridine is a anhydrous nucleoside; 2'-O-Methyl nucleoside.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.163 mL	20.8151 mL	41.6302 mL
5 mM	0.8326 mL	4.163 mL	8.326 mL
10 mM	0.4163 mL	2.0815 mL	4.163 mL
50 mM	0.0833 mL	0.4163 mL	0.8326 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481