


16,17-Epoxypregnenolone acetate

Chemical Properties

CAS No. :	34209-81-9
Formula:	C ₂₃ H ₃₂ O ₄
Molecular Weight:	372.5
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	16,17-Epoxypregnenolone acetate (16alpha,17alpha-Epoxypregnenolone acetate) can be produced by Burkholderia cepacia SE-1.
Targets(IC50)	Others

Solubility Information

Solubility	DMSO: 6 mg/mL (16.11 mM), Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.6846 mL	13.4228 mL	26.8456 mL
5 mM	0.5369 mL	2.6846 mL	5.3691 mL
10 mM	0.2685 mL	1.3423 mL	2.6846 mL
50 mM	0.0537 mL	0.2685 mL	0.5369 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Zhu X, et al. Biomed Res Int. 2016;2016:5727631.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481