


## Sulfo-SIAB sodium

## Chemical Properties

CAS No. :	144650-93-1
Formula:	C13H10IN2NaO8S
Molecular Weight:	504.19
Storage:	Keep away from direct sunlight Powder: -20°C for 3 years   In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



## Biological Description

Description	Sulfo-SIAB sodium is a non-cleavable monovalent bifunctional linker.
Targets(IC50)	ADC Linker,PROTAC Linker

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.9834 mL	9.9169 mL	19.8338 mL
5 mM	0.3967 mL	1.9834 mL	3.9668 mL
10 mM	0.1983 mL	0.9917 mL	1.9834 mL
50 mM	0.0397 mL	0.1983 mL	0.3967 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

## Reference

Gubala V, et al. A comparison of mono and multivalent linkers and their effect on the colloidal stability of nanoparticle and immunoassays performance. Talanta. 2010 Jun 15;81(4-5):1833-9.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481