

Isouvaretin

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Isouvaretin demonstrates potent antibacterial activity against B. subtilis and S. epidermidis.
Targets(IC50)	Others,Antibacterial

Reference

Maeda G, et al. Oxygenated Cyclohexene Derivatives from the Stem and Root Barks of Uvaria pandensis. J Nat Prod. 2021 Nov 21.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481