

Compound T129237

Chemical Properties


CAS No. :

Formula: C₃₄H₃₆O₁₅

Molecular Weight: 684.65

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	(1S,2S,4S,5R,6R,7S,9R,12R)-12-(acetyloxy)-6-[(acetyloxy)methyl]-5,7-bis(furan-3-carboxyloxy)-2-hydroxy-2,10,10-trimethyl-11-oxatricyclo[7.2.1.0 ^{1,6}]dodecan-4-yl furan-3-carboxylate is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.4606 mL	7.303 mL	14.606 mL
5 mM	0.2921 mL	1.4606 mL	2.9212 mL
10 mM	0.1461 mL	0.7303 mL	1.4606 mL
50 mM	0.0292 mL	0.1461 mL	0.2921 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481