

## Oxaprozin D4

## Chemical Properties


CAS No. :

Formula: C<sub>18</sub>H<sub>11</sub>D<sub>4</sub>NO<sub>3</sub>

Molecular Weight: 297.34

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description

Oxaprozin-D4 is the deuterium labeled Oxaprozin (TMSM-1824), is a non-steroidal anti-inflammatory drug (NSAID).

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.3632 mL	16.8158 mL	33.6315 mL
5 mM	0.6726 mL	3.3632 mL	6.7263 mL
10 mM	0.3363 mL	1.6816 mL	3.3632 mL
50 mM	0.0673 mL	0.3363 mL	0.6726 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481