

Propargyl-PEG2-CH₂COOH

Chemical Properties

CAS No. : 944561-46-0


Formula: C₉H₁₄O₅

Molecular Weight: 202.2

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	Propargyl-PEG2-CH ₂ COOH is a PEG-based linker for PROTACs (joining two essential ligands) that facilitates selective protein degradation by leveraging the ubiquitin-proteasome system within cells.
Targets(IC50)	Others,PROTAC Linker

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.9456 mL	24.728 mL	49.456 mL
5 mM	0.9891 mL	4.9456 mL	9.8912 mL
10 mM	0.4946 mL	2.4728 mL	4.9456 mL
50 mM	0.0989 mL	0.4946 mL	0.9891 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Guangrong Zheng, et al. Compounds that induce degradation of anti-apoptotic bcl-2 family proteins and the uses thereof. WO2017184995A1.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481