

Resorufin pentyl ether

Chemical Properties

CAS No. : 87687-03-4


Formula: C₁₇H₁₇NO₃

Molecular Weight: 283.32

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	fluorometric cytochrome P450 substrate
Targets(IC50)	Others,Antibacterial

Solubility Information

Solubility	DMSO: <1 mg/mL, Ethanol: <1 mg/mL, (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	----------------------------------------------------------------------------------------------------------

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.5296 mL	17.6479 mL	35.2958 mL
5 mM	0.7059 mL	3.5296 mL	7.0592 mL
10 mM	0.353 mL	1.7648 mL	3.5296 mL
50 mM	0.0706 mL	0.353 mL	0.7059 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481