


## Erythartine

## Chemical Properties

CAS No. :	51666-26-3
Formula:	C <sub>19</sub> H <sub>23</sub> NO <sub>4</sub>
Molecular Weight:	329.39
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Erythartine is a natural product of Erythrina, Fabaceae. The catalog number is TN3996 and the CAS number is 51666-26-3. Erythartine can be used as a reference standard.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.0359 mL	15.1796 mL	30.3591 mL
5 mM	0.6072 mL	3.0359 mL	6.0718 mL
10 mM	0.3036 mL	1.518 mL	3.0359 mL
50 mM	0.0607 mL	0.3036 mL	0.6072 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481