


3'',4''-Di-O-p-coumaroylquercitrin

Chemical Properties

| | |
|-------------------|---|
| CAS No. : | 437615-43-5 |
| Formula: | C39H32O15 |
| Molecular Weight: | 740.67 |
| Storage: | Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA. |



Biological Description

| | |
|-------------|--|
| Description | 3'',4''-Di-O-p-coumaroylquercitrin is a natural product for research related to life sciences. The catalog number is TN5487 and the CAS number is 437615-43-5. |
|-------------|--|

Preparing Stock Solutions

| | 1mg | 5mg | 10mg |
|-------|-----------|-----------|------------|
| 1 mM | 1.3501 mL | 6.7506 mL | 13.5013 mL |
| 5 mM | 0.270 mL | 1.3501 mL | 2.7003 mL |
| 10 mM | 0.135 mL | 0.6751 mL | 1.3501 mL |
| 50 mM | 0.027 mL | 0.135 mL | 0.270 mL |

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481