


## Secologanin dimethyl acetal

## Chemical Properties

CAS No. :	77988-07-9
Formula:	C <sub>19</sub> H <sub>30</sub> O <sub>11</sub>
Molecular Weight:	434.43
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Secologanin dimethyl acetal induces significant neurite outgrowth.
Targets(IC50)	Others,Antibacterial

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.3019 mL	11.5093 mL	23.0187 mL
5 mM	0.4604 mL	2.3019 mL	4.6037 mL
10 mM	0.2302 mL	1.1509 mL	2.3019 mL
50 mM	0.046 mL	0.2302 mL	0.4604 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481