


## 3'-Deoxy-N6-methyladenosine

## Chemical Properties

CAS No. :	3616-27-1
Formula:	C <sub>11</sub> H <sub>15</sub> N <sub>5</sub> O <sub>3</sub>
Molecular Weight:	265.27
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	3'-Deoxy-N6-methyladenosine is a Nucleoside Derivative - 3'-Modified nucleoside, 6-Modified purine nucleoside;3'-Deoxy nucleoside.
Targets(IC50)	Nucleoside Antimetabolite/Analog,Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.7697 mL	18.8487 mL	37.6974 mL
5 mM	0.7539 mL	3.7697 mL	7.5395 mL
10 mM	0.377 mL	1.8849 mL	3.7697 mL
50 mM	0.0754 mL	0.377 mL	0.7539 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481