


## 3-Epiturraeanthin

## Chemical Properties

CAS No. : 22415-24-3  
 Formula: C<sub>32</sub>H<sub>50</sub>O<sub>5</sub>  
 Molecular Weight: 514.747  
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
 Actual storage temperature shall be subject to the COA.



## Biological Description

Description	3-Epiturraeanthin is a natural product for research related to life sciences. The catalog number is TN2932 and the CAS number is 22415-24-3.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.9427 mL	9.7135 mL	19.4269 mL
5 mM	0.3885 mL	1.9427 mL	3.8854 mL
10 mM	0.1943 mL	0.9713 mL	1.9427 mL
50 mM	0.0389 mL	0.1943 mL	0.3885 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481