


TCS 2312

Chemical Properties

CAS No. : 838823-31-7
 Formula: C₂₅H₂₄N₄O₂
 Molecular Weight: 412.48
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	checkpoint kinase 1 (chk1) inhibitor
Targets(IC50)	Others,Chk

Solubility Information

Solubility	DMSO: <8.25 mg/mL, Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.4244 mL	12.1218 mL	24.2436 mL
5 mM	0.4849 mL	2.4244 mL	4.8487 mL
10 mM	0.2424 mL	1.2122 mL	2.4244 mL
50 mM	0.0485 mL	0.2424 mL	0.4849 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481