

Compound T130913

Chemical Properties


CAS No. :

Formula: C₂₂H₃₆O₄

Molecular Weight: 364.53

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	5-[(1S,2R,4aR)-5-[(acetyloxy)methyl]-1,2,3,4,4a,7,8,8a-octahydronaphthalen-1-yl]-3-methylpentanoic acid is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.7433 mL	13.7163 mL	27.4326 mL
5 mM	0.5487 mL	2.7433 mL	5.4865 mL
10 mM	0.2743 mL	1.3716 mL	2.7433 mL
50 mM	0.0549 mL	0.2743 mL	0.5487 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481